

Tutorial of KGMN result visualization and analysis

Zhiwei Zhou

2022-06-05

Introduction

Unknown metabolite annotation is one of long-standing challenges in untargeted metabolomics. We developed an approach, namely, knowledge-guided multi-layer network (KGMN), to enable global metabolite annotation from knowns to unknowns in untargeted metabolomics. The KGMN approach integrated three-layer networks, including knowledge-based metabolic reaction network (Network 1), knowledge-guided MS/MS similarity network (Network 2), and global peak correlation network (Network 3). This tutorial will help users to visualize, reproduce and investigate putatively annotated known and unknown metabolites from KGMN.

1. Installation

The analysis and visualization of KGMN results mainly relies on R package – MetDNA2Vis, and its depended R packages; The Cytoscape software is used for manually visualize networks, and interactively investigate results of KGMN; The ChemDraw software is involved for drawing chemical structures.

- Install R packages

```
# Install related packages
if(!require(devtools)){
install.packages("devtools")
}

if(!require(BiocManager)){
install.packages("BiocManager")
}

# Install CRAN/Bioconductor packages
required_pkgs <- c("dplyr","tidyr","readr","CHNOSZ","igraph",
  "magrittr","ggplot2","ggraph","tidygraph")
list_installed <- installed.packages()

new_pkgs <- required_pkgs[!(required_pkgs %in% list_installed[, 'Package'])]
```

```
if (length(new_pkgs) > 0) {
  BiocManager::install(new_pkgs)
} else {
  cat('Required CRAN/Bioconductor packages installed\n')
}
```

Install ZhuLab packages

```
devtools::install_github("ZhuMetLab/SpectraTools")
```

```
devtools::install_github("ZhuMetLab/MetDNA2Vis")
```

- Cytoscape software (Version 3.8 or higher required): <https://cytoscape.org/>
- ChemDraw software (Version 19.0 or higher required): <https://perkinelmerinformatics.com/products/research/chemdraw>

2. Step-by-step instruction for visualization

In this part, we introduce how to visualize multi-layer networks from KGMN. It will help users to reproduce figures in KGMN manuscripts. Here, the Human NIST urine (Positive data, used in KGMN manuscript) was used as a demo dataset. This data set have been processed and exported by **MetDNA2 web server** (version 1.0.4). The raw data files and results can be downloaded at [here](#). The more details of sample extraction and data preprocessing can be found in our KGMN manuscript.

The step-by-step demonstration is provided as below.

2.1 Download demo data and unzip the archive.

- All required intermediate files for visualization is provided in '06_visualization' folder.

名称	修改日期	类型	大小
00_annotation_table	2022/6/4 15:36	文件夹	
02_result_MRN_annotation	2022/6/4 15:36	文件夹	
04_biology_intepretation	2022/6/4 15:36	文件夹	
05_analysis_report	2022/6/4 15:36	文件夹	
06_visualization	2022/6/4 15:30	文件夹	
data.csv	2022/1/17 9:12	Microsoft Excel ...	
NIST_urine01_pos-NIST_urine01.mgf	2022/1/17 9:10	MGF 文件	
NIST_urine02_pos-NIST_urine02.mgf	2022/1/17 9:12	MGF 文件	
NIST_urine03_pos-NIST_urine03.mgf	2022/1/17 9:12	MGF 文件	
NIST_urine04_pos-NIST_urine04.mgf	2022/1/17 9:10	MGF 文件	
para_list.txt	2022/6/4 15:33	文本文档	2 KB
QC_pos-QC.mgf	2022/1/17 9:12	MGF 文件	9,687 KB
RT_recibration_table.csv	2022/1/17 9:12	Microsoft Excel ...	1 KB
sample.info.csv	2022/1/17 9:12	Microsoft Excel ...	1 KB

list_peak_group	2022/4/13 3:47	文件	151 KB
list_peak_group_annotation_concis...	2022/4/13 3:28	R Workspace	6,730 KB
ms2_data.msp	2022/4/13 3:20	Windows Install...	1,318 KB
ms2_data.RData	2022/4/13 3:21	R Workspace	591 KB
peak_group_id_table	2022/4/13 3:45	文件	38 KB
table_identification	2022/4/13 3:47	文件	138 KB

2.2 Preparing.

- Set the working directory ('your_path/06_visualization') and load required packages. Then, please check required files whether existed.

```
# load packages
library(MetDNA2Vis)
library(CHNOSZ)
library(dplyr)

# check required files
checkFiles4Vis()

## Check required files ...
## Check required files: done!
```

2.3 Reconstruct and export global multi-layer networks.

2.3.1 Network 1

The network 1 is the knowledge-guided metabolic reaction network. For knowns, the KEGG reaction pair network is directly used. For unknowns, an extended KEGG reaction pair network was used. The network expansion is performed with in-silico enzymic reactions (via Biotransformer), and further connected with KEGG reaction pair network. The details of network construction and expansion are described in our KGMMN manuscript. It should be noted that the KEGG reaction pair network and extended network were built in advance.

To export the network 1, it is easy to run reconstructNetwork1 function as below:

```
# export network 1 for visualization
reconstructNetwork1(is_unknown_annotation = TRUE)
```

The network files will be exported in '00_network1' folder. It contains two files, including "edge_table.tsv" and "node_table.tsv" (Figure 2.3.1). These tables can be imported into Cytoscape software for visualization.

名称	修改日期	类型	大小
00_files_network1			
01_files_network2			
02_files_network3			
03_subnetworks			
.Rhistory	2022/6/5 10:51	文件夹	
list_peak_group	2022/6/5 12:12	R History 源文件	2 KB
list_peak_group_annotation_concis...	2022/4/13 3:47	文件	151 KB
ms2_data.msp	2022/4/13 3:28	R Workspace	6,730 KB
ms2_data.RData	2022/4/13 3:20	Windows Install...	1,318 KB
peak_group_id_table	2022/4/13 3:21	R Workspace	591 KB
table_identification	2022/4/13 3:45	文件	38 KB
	2022/4/13 3:47	文件	138 KB

2.3.2 Network 2


The network 2 is a knowledge-guided MS/MS network. Although it calls MS/MS network, differing to MS/MS network (mainly based on MS2), the linkage (edge) of network2 has a prerequisite. It requires a reasonable reaction relationship and definitive structure candidate first. As a result, their retention time can also be predicted. In other words, two linked nodes indicate 4 messages. Their candidates of these nodes have (1) reasonable reaction relationships, (2) low m/z errors, (3) low RT error against with predicted RT values, and (4) MS/MS similarity. It should be note that optimized network2 required to be reconstructed from KGMN exported results, because the global peak correlation network remove and collapse some error nodes and edges in prior analysis. This process usually requires 10-20 min to complete.

To export the network 2, it is easily to run reconstructNetwork2 function as below:

```
# Modify format of KGMN result
annotation_table <- reformatTable1()

# Export global network2 files
reconstructNetwork2(annotation_table = annotation_table,
  is_unknown_annotation = TRUE)
```

The networks files will be exported in '01_network2' folder. The "edge_table.tsv" and "node_table.tsv" in this folder can be imported to Cytoscape.



名称	修改日期	类型	大小
edge_table_network2.tsv	2022/6/4 21:22	TSV 文件	165 KB
network2_obj.RData	2022/6/4 21:22	R Workspace	68 KB
network2_result.RData	2022/6/4 21:22	R Workspace	59 KB
node_table_network2.tsv	2022/6/4 21:22	TSV 文件	80 KB

名称	修改日期	类型	大小
.Rhistry	2022/6/5 12:12	R History 源文件	2 KB
list_peak_group	2022/4/13 3:47	文件	151 KB
list_peak_group_annotation_concis...	2022/4/13 3:28	R Workspace	6,730 KB
ms2_data.msp	2022/4/13 3:20	Windows Install...	1,318 KB
ms2_data.RData	2022/4/13 3:21	R Workspace	591 KB
peak_group_id_table	2022/4/13 3:45	文件	38 KB
table_identification	2022/4/13 3:47	文件	138 KB


2.3.3 Network 3

The network 3 is the global peak correlation network. This network recognized abiotic peaks derived from peaks from network 2, including adducts, isotopes, neutral losses, and in-source fragments (ISF). The network 3 is used to optimize the annotation and linkage of network 2. The optimization has been completed in KGMN analysis. The details of network 3 construction and optimization can be found in our manuscript.

To export the network 3, it is easily to run reconstructNetwork3 function as below:

```
# export network3
reconstructNetwork3()
```

The networks files will be exported in '**02_files_network3**' folder. The "edge_table.tsv" and "node_table.tsv" in this folder can be imported to Cytoscape for visualization.



名称	修改日期	类型	大小
edge_table_network3.tsv	2022/6/4 21:41	TSV 文件	294 KB
network3_obj.RData	2022/6/4 21:41	R Workspace	151 KB
network3_result.RData	2022/6/4 21:41	R Workspace	113 KB
node_table_network3.tsv	2022/6/4 21:41	TSV 文件	340 KB


名称	修改日期	类型	大小
.Rhistry	2022/6/5 12:12	R History 源文件	2 KB
list_peak_group	2022/4/13 3:47	文件	151 KB
list_peak_group_annotation_concis...	2022/4/13 3:28	R Workspace	6,730 KB
ms2_data.msp	2022/4/13 3:20	Windows Install...	1,318 KB
ms2_data.RData	2022/4/13 3:21	R Workspace	591 KB
peak_group_id_table	2022/4/13 3:45	文件	38 KB
table_identification	2022/4/13 3:47	文件	138 KB

2.4 Visualize global networks with Cytoscape


Above networks (Network 1-3) can be imported to Cytoscape software tool for visualization. The process of network visualization is generally similar. Here, we use the above network 1 as a demonstration. The version of Cytoscape used here is 3.8.2.

Below is the step-by-step instruction:


1. **Import edge file.** Select the “edge_table.tsv” file and open it in the box.




2. **Assign column attributes.** Click the “from” column and select it as “source node”. Similarly, click the “to” column and select it as “target node”. After assigning attributes, click **OK** to construct a network.




3. **Import node file.** Select the “node_table.tsv” file and open it in the box.




4. Select the “name” column as a key. Then, click the **OK** button.



5. **Modify the style for visualization.** Click the Style type, you can adjust node shapes and colors, edge types and colors.




To help users reproduce our plot quickly, users can directly import our style file. The styles of different networks are provided [here](#).



2.5 Select and export interesting subnetwork

Through above procedures, users can easily visualize global network 1-3. With such global networks, users can find interesting subnetworks in Cytoscape. The Cytoscape supports a interactively investigation. **It should be note that the targeted subnetwork selection is customized.** Users can directly find interesting nodes from KGMN annotation results, or considering more information, like in-silico MS/MS, chemical structure and/or statistics analysis. For example, in KGMN manuscript, we combined MASST to select an unknown subnetwork of M262T526 (**Figure 5e in manuscript**). This unknown peak is putatively annotated as O-sulfotyrosine, and this annotation is from M182T541-Tyrosine. This subnetwork consists of 2 peaks and 2 metabolites. Here, we mainly introduce how to export and visualize this subnetwork. First, export network 1 of this subnetwork. **Note:** the export and visualization require intermediate results from global networks. Therefore, please run global peaks export first. To export the subnetwork 1, please directly run retrieveSubNetwork1 function as below.

```
# network 1 of unknown peak subnetwork
# Note: the folder_output should keep same among different layer subnetworks
retrieveSubNetwork1(centric_met = c('C00082', 'KeggExd000923'),
  is_unknown_annotation = TRUE,
  folder_output = c('M182T541_M262T526'))
```




The networks files will be exported in '03_subnetworks/your_defined_folder/network 1' folder. Here, the exported folder is "M182T541_M262T526". The "edge_table.tsv" and "node_table.tsv" in this folder can be imported to Cytoscape for visualization. **Note:** if you run in RStudio, the preview plot of subnetwork 1 will be directly shown in the plot panel.

Similarly, export network 2 and network 3 of this subnetwork can be completed through running `retrieveSubNetwork2` and `retrieveSubNetwork3` functions, respectively. The preview plots of subnetwork 2 and subnetwork 3 will be shown in the plot panel if you run in RStudio.


```
# network 2 of unknown peak subnetwork
retrieveSubNetwork2(from_peak = 'M182T541',
  end_peak = 'M262T526',
  folder_output = c('M182T541_M262T526'))

## Using `sugiyama` as default layout
```



```
# network 3 of unknown peak subnetwork
retrieveSubNetwork3(base_peaks = c('M182T541', 'M262T526'),
  base_adducts = c('[M+H]+', '[M+H]+'),
  folder_output = c('M182T541_M262T526'))


## Using `stress` as default layout
```



The network 2 and network 3 of the subnetwork can be further merged through running mergeSubnetwork function. The 'network_merge' folder contains node table and edge table for reproduce the merged network.

```
# merge subnetwork
mergeSubnetwork(from_peak = 'M182T541',
  end_peak = 'M262T526',
  folder_output = 'M182T541_M262T526')

## Using `stress` as default layout
```



Finally, the folder of subnetwork is organized like below. Each folder contains related files of each network for further visualization in other tools (e.g. Cytoscape).

00_files_network1	2022/6/5 10:47	文件夹	
01_files_network2	2022/6/4 21:22	文件夹	
02_files_network3	2022/6/4 21:41	文件夹	
03_subnetworks	2022/6/5 10:51	文件夹	
.Rhistory	2022/6/5 12:12	R History 源文件	2 KB
list_peak_c	名称	修改日期	类型 大小
list_peak_c	M182T541_M262T526	2022/6/5 11:44	文件夹 6,730 KB
ms2_data.msp	2022/4/13 3:20	Windows Install...	1,318 KB
ms2_data.RData	2022/4/13 3:21	R Workspace	591 KB
peak_grou	名称	修改日期	类型 大小
table_ider	network_merge	2022/6/5 11:44	文件夹
	network1	2022/6/5 10:51	文件夹
	network2	2022/6/5 11:27	文件夹
	network3	2022/6/5 11:30	文件夹

3. The script for visualization

Here is a script contains above codes to help to reproduce above analysis quickly.

```
# load packages
library(CHNOSZ)
library(dplyr)
```

```

library(MetDNA2Vis)

# set working directory
setwd('D:/project/00_zhulab/01_metdna2/00_data/20220602_visualization_kgmn/Demo_MetDNA2_NI
ST_urine_pos/06_visualization/')

# Export global networks
# construct network 1
reconstructNetwork1(is_unknown_annotation = TRUE)

# construct network 2
annotation_table <- reformatTable1()
reconstructNetwork2(annotation_table = annotation_table)

# construct network 3
reconstructNetwork3()

# Export subnetworks -----
# network 1 of unknown peak subnetwork
# Note: the folder_output should keep same among different layer subnetworks
retrieveSubNetwork1(centric_met = c('C00082', 'KeggExd000923'),
  is_unknown_annotation = TRUE,
  folder_output = c('M182T541_M262T526'))

# network 2 of unknown peak subnetwork
retrieveSubNetwork2(from_peak = 'M182T541',
  end_peak = 'M262T526',
  folder_output = c('M182T541_M262T526'))

# network 3 of unknown peak subnetwork
retrieveSubNetwork3(base_peaks = c('M182T541', 'M262T526'),
  base_adducts = c('[M+H]+', '[M+H]+'),
  folder_output = c('M182T541_M262T526'))

# merge subnetwork
mergeSubnetwork(from_peak = 'M182T541',
  end_peak = 'M262T526',
  folder_output = 'M182T541_M262T526')

```